CBC NMR Laboratory Reopening

April 19, 2021

- The maximum number of occupancies is 3 and 6 for 045 BRL and 049 BRL, respectively.
- 045 Brown Lab is the designated entrance to the NMR Facility
- You must wash your hands before entering the NMR Facility at the sink outside 045 Brown Lab



- Follow the flow of traffic designated by the red arrows on the floor.
- You must sanitize sample depth-gage, auto-sample touch pad, mouse, keyboard, and desktop with 70% IPA on a paper towel after submitting your samples
- The door by the AVIII 400 (manual) is the designated exit







Please also note:

- The NMR Lab will be closed from 8:30 to 11:30 on Monday mornings for cryogen fills and lab disinfection.
- The NMR Lab or a part of the Lab will be closed during training sessions (TBA)